

GIENIC | PureWick™

Urine Collection System

An innovative way to manage urinary incontinence by using gentle suction and an external catheter to draw urine into a sealed canister.



Get a better night's sleep.



Keep skin dry + reduce incontinence-associated dermatitis.



Less high-risk, night-time bathroom trips.

Male External Catheter

PureWick™

A non-invasive option for urine output management for men.

When the user urinates on the external catheter, the PureWick™ Urine Collection system uses suction to pull the voided urine through tubing to the connected collection canister, where it is stored until ready for proper disposal.



PureWick™ External Catheters are connected to the PureWick™ Urine Collection system. They sit outside the body, where they are held in place by an adhesive seal for men.

“By using the PureWick™ System, my patients have experienced better sleep. It also allows their caregivers and family members to have better sleep since they are not having to get up multiple times overnight to assist their loved ones with going to the restroom and with changing wet incontinence garments.

- Dr. Amanda Lathia; MD, Geriatric Medicine



Female External Catheter

PureWick™

A non-invasive, dignified urine output management option for female.

PureWick™ Female External Catheter is designed to fit around the contours of your body for a comfortable fit that stays in place.



A breakthrough external catheter for managing urinary incontinence at home.




No bathroom trips during the night, and far better sleep - what a relief!

- Maryann B.

GIENIC

Enquire today for your complete Aged Care + Clinical Supplies.



Visit www.gienic.com.au/purewick for more details.

1300 00 22 90 | sales@gienic.com.au